

Cerenkov Imaging of Radioisotopes in IVIS systems

Author


Chaincy Kuo, Ph.D.
Staff Scientist, In Vivo Imaging
PerkinElmer, Inc.
68 Elm Street
Hopkinton, MA 01748



Cerenkov Emission from radioisotopes in tissue

Optical imaging detects photons in the visible range of the electromagnetic spectrum. PET and SPECT imaging instruments are sensitive to photons in the much higher energy range of x-rays and gamma rays. While the PET and SPECT probes which can generate Cerenkov light in tissue will continue to produce the relevant gamma- and x-rays, visible photons will be produced from the Cerenkov emission, which the IVIS® will detect.

In beta decay emitters such as PET probes and isotopes such as ^{90}Y , ^{177}Lu , ^{131}I and ^{32}P , the beta particle will travel in the tissue until it either annihilates with an electron or loses momentum due to viscous electromagnetic forces.




It is possible that the beta (electron or positron) is relativistic, traveling faster than the speed of light in the tissue. While it is impossible to travel than the speed of light in a vacuum (c), the speed of light in tissue is $v = c / n$, where n is the tissue index of refraction and $n \geq 1$. Cerenkov photons will be generated by a relativistic charged particle in a dielectric medium such as tissue.

Highly Sensitive CCD's

As these photons are in the visible wavelength band, the imaged photons are subject to scattering and absorption by the tissue to the degree that their propagation can be approximated as a diffuse process. The photons eventually reach the animal surface and are imaged by the highly sensitive CCD.

IVIS CCD's are optimal for imaging the low photon counts generated by Cerenkov emission from radioisotopes as the CCD's are cooled to -90°C , and back-thinned with a high quantum efficiency ($>80\%$) over $450\text{ nm} - 750\text{ nm}$. These characteristics offer a low read noise of $\sim 3\text{ e RMS}$ and a dark current of $< 4 \times 10^{-4}\text{ e/sec/pixel}$ at -90°C . The high resolution CCD's are 1024 or 2048 pixels on a side with $13\text{ }\mu\text{m}$ pixel size. For improved sensitivity necessary for Cerenkov imaging, the pixels can be binned up to 16×16 , and integration times range from 0.5 seconds to 10 minutes.



Radioisotope biodistribution optical imaging with DyCE


Biodistribution curves of radioisotopes can be determined by imaging in the IVIS and analyzing the images in DyCE™ software.

^{18}F -FDG kinetics

A mouse bearing a subcutaneous 4T1-luc2 tumor in its right flank was injected with $315\text{ }\mu\text{Ci}$ of ^{18}F -FDG intravenously. 55 seconds post-injection, the animal was imaged dynamically starting 55 seconds post-injection to capture the distribution of ^{18}F -FDG in the mouse body via Cerenkov light from positron emission. No luciferin was administered to the mouse to avoid contamination of luciferase signal into the Cerenkov signal.

DyCE was applied to the series of 2D time-domain 'Open' filter images of Cerenkov light emission from positron emission. DyCE characterizes the temporal change in pixel intensities of the 2D images.


Without a priori identification of organ regions, the unmixed component amplitudes map to pixel regions in the 2D image which can be associated with anatomical features. By corresponding pixel color-mapping to component temporal evolution, kinetics of ^{18}F -FDG can be characterized.



⁹⁰Y-AABD kinetics

Mice injected with 10 μ Ci of ⁹⁰Y-AABD and imaged dynamically in optical imaging system starting 55 seconds post-injection. Imaging parameters: 60 second exposures at bin level 16 and f-stop 1, dorsal imaging using the 'Open' filter.

⁹⁰Y is a beta-minus emitting isotope where optical imaging of Cerenkov radiation is the only means for in vivo characterization of biodistribution in small animals. DyCE was performed to determine kinetics of ⁹⁰Y-AABD. Relative abundances of tracer uptake are associated in the kidneys, muscle and brain regions in the 2D component overlay image.



Summary

Cerenkov imaging of radiopharmaceuticals injected in small animals is a low-cost solution to PET and SPECT instrumentation.

IVIS systems have highly sensitive CCD's required for the low signal Cerenkov light, and offer high signal to noise ratios.

High animal throughput is also available in all IVIS systems, where up to 5 animals can be imaged in 2D simultaneously.

Kinetics of PET probes and radio-therapeutics such as ⁹⁰Y and ¹³¹I can be followed with optical imaging and analyzed with DyCE.

Table 1. List of radioisotopes produced and distributed by PerkinElmer. Those which will produce Cerenkov Light in vivo are highlighted in green. Through its worldwide distribution capabilities, PerkinElmer can assure delivery of materials to many areas of the world soon after production. If you require additional information on a specific isotope or information on a product that is not listed below, please contact us at Radiotherapeutics@perkinelmer.com.

Isotope	Half-life	Radiation	Kinetic Energy (keV)	Fresh Lot Date	Application
Indium-111 No Carrier Added Indium chloride in 0.05M HCl	2.83 days	Gamma (γ)	171 245	Every Friday	In Vitro assays In Vivo Imaging-SPECT
Iodine-124 No Carrier Added Sodium Iodide in 20 mM NaOH pH >10	4.2 days	Beta+ (β^+)	686 974	Every Monday	In Vivo Imaging-Cerenkov Immuno-PET
Iodine-125 No Carrier Added Sodium Iodide in NaOH pH 12-14 or 8-11	59.4 days	Gamma (γ)	35.50	Every Wednesday	In Vitro assays In Vivo Imaging-SPECT
Iodine-131 No Carrier Added Sodium Iodide in NaOH pH 12-14	8.04 days	Gamma (γ) Beta- (β^-)	364 637 284 606	Every Monday	β^- Therapy In Vivo Imaging-Cerenkov & SPECT
Phosphorus-32 No Carrier Added Please Inquire	14.29 days	Beta- (β^-)	1710	Every Wednesday	In Vitro assays In Vivo Imaging-Cerenkov
Yttrium-90 No Carrier Added Yttrium chloride in 0.05M HCl >500mCi/ml	64.1 hours	Beta- (β^-) Beta+ (β^+)	2.28	Multiple times per week	β^- Therapy PET Calibration In Vivo Imaging-Cerenkov & PET
Zirconium-89 No Carrier Added in 1M Oxalic Acid pH <4	78.41 hours	Beta+ (β^+)	396	Every Thursday	In Vivo Imaging-Cerenkov & Immuno-PET
Lutetium-177 ~27Ci/mg Lutetium Chloride In 0.05M HCl ~3Ci/ml	6.64 days	Beta- (β^-) Gamma (γ)	0.50 208	Every Week	β^- Therapy In Vivo Imaging-SPECT Cerenkov

Selected References

Jelley, J (1958). *Cerenkov Radiation and its Applications*. New York, NY: Pergamon Press, Inc.

Robertson R, Germanos MS, Li C, Mitchell GS, Cherry SR and Silva MD. (2009) Optical imaging of Cerenkov light generation from positron emitting radiotracers. *Phys Med Biol*, 54:N355-365.

Cheng Z, Liu HG, Ren G, Liu SL, Zhang XF, Chen LX and Han PZ. (2010) Optical imaging of reporter gene expression using a positron-emission tomography probe. *J Biomed Opt*, 15:060505.

Li, C, Mitchell, GS, and Cherry, SR (2010). Cerenkov luminescence tomography for small-animal imaging. *Optics Letters*, 35 (7), 1109-1111.

Liu H, Ren G, Miao Z, Zhang X, Tang X, Han P, Gambhir SS and Cheng Z. Molecular Optical Imaging with Radioactive Probes. *PLoS One* 2010; 5: e9470.

Ruggiero A, Holland JP, Lewis JS and Grimm J. (2010) Cerenkov luminescence imaging of medical isotopes. *J Nucl Med*, 51: 1123-1130.

Kuo C, Mitchell GS, Gill RK, Boucher DL, Rice BW, and Cherry SR (2011), Multivariate curve estimation of radiotracer kinetics by Cerenkov radiation imaging in small animals, poster presentation at *World Molecular Imaging Congress 2011*, P065.

Spinelli AE, Kuo C, Rice BW, Calandrino R, Marzola P, and Sbarbati A (2011). Multispectral Cerenkov luminescence tomography for small animal optical imaging. *Optics Express*, 12605-12618.

Mitchell GS, Gill RK, Boucher DL, Li C, and Cherry SR (2011). In vivo Cerenkov luminescence imaging: a new tool for molecular imaging. *Philosophical Transactions of the Royal Society A*, 4605-4619.

Kuo C, Mitchell GS, Gill RK, Rice BW, Cherry SR (2012), Radionuclide quantification with multi-spectral Cerenkov Tomography for a small animal imaging system, poster presentation at *World Molecular Imaging Congress 2012*, P380.

Zhang X, Kuo C, Moore A, Ran C, (2013), In Vivo Optical Imaging of Interscapular Brown Adipose Tissue with (18)F-FDG via Cerenkov Luminescence Imaging, *PLoS One*, 8:e62007.