

UWCCC Biotech Attune Nxt V6 Blank Instrument Map River

488 Blue 2 Detectors		633 Red 3 Detectors			405 Violet 6 Detectors						561 YellowGreen 3 Detectors			
Detector	BL1	BL2	RL1	RL2	RL3	VL1	VL2	VL3	VL4	VL5	VL6	YL1	YL2	YL3
Filter	530/30	695/40	670/14	720/30	780/60	450/40	525/50	610/20	660/20	710/50	780/60	585/16	620/15	780/60
Dichroic	503LP	555LP	654LP	690LP	740LP	417LP	495LP	555LP	635LP	680LP	740LP	577LP	600LP	650LP

Notes:

Side Scatter(SSC) Off Violet Possible With No Lyse/No Wash Filter Set	Fluorochromes are listed for optimal excitation and emission.
GFP/YFP/mCherry; Fluorecent Protein Filter Set Is Available For Optimal Configuration	- BV, BUV, BB dyes require the use of the Brilliant Dye Staining Buffer.
PerCP alone does not work on the Attune (goes into "triplet state" and stops fluorescing); PerCP tandems are okay.	- SuperBright dyes require the use of the SuperBright Staining Buffer.
*Starred fluorochromes will show up in all channels listed, but are brightest in the Bolded channel.	- NovaFluor dyes require the use of the Cellblox Plus Blocking Buffer.
Red Text Indicates Filter Swap. Optimal Filter Listed. Confirm Filter Availability with UWCCC Flow Lab Staff.	- NovaFluor dyes can NOT be used with DNA dyes (DAPI, 7AAD, PI, etc.).
For newly released or less common fluorochromes please see spectra viewers for usability on this instrument or email uwflow@uwcarbone.wisc.edu for guidance.	